

## LASMIK-VET Kit for Veterinary Medicine



Initially, the kit was developed to specifically treat and rehabilitate horses and maintain their fitness. Still, the capabilities of the devices are much greater, allowing for the most effective treatment of other domestic and farm animals.

The kit includes two low-level laser therapy devices with an autonomous battery power supply

### «LASMIK-VET»

The device is a matrix pulsed laser, an analogue of the ML-904-80 laser emitting head for Matrix and Lasmik laser therapeutic devices. It is intended for local exposure.

- wavelength – 904 nm;
- operating mode – pulsed;
- light pulse duration – 100 ns;
- power – 80 W;
- number of laser diodes – 8 pcs.;
- frequencies: 80 and 10,000 Hz;
- timer: 0.5, 1, 2, 5 min.

### «LASMIK-AP»

It is designed for laser acupuncture and local illumination.

Main technical characteristics of the device:

- wavelength – 635 nm;
- power – 5–10 mW;
- number of laser diodes – 1 pc.;
- mode of operation – continuous work;
- timer – 30 s.

### Optical attachments

- Acupuncture nozzle A-3 – 1 pc.
- PMN nozzle (transparent protective) – 1 pc.
- Mirror attachment ZN-35 – 1 pc.
- MM-50 attachment (magnetic 50 mT) – 1 pc.
- ZM-50 attachment (magnetic 50 mT) – 1 pc.

### Additionally

- Charger – 1 pc.
- Charging cord – 2 pcs.
- Book “Low-level laser therapy for sports horses” (in Russian and English) – 1 set.
- Book “Low-level laser therapy for pets” (in Russian and English) – 1 set.
- Bag – 1 pc.

Quite often, customers also purchase LASMIK<sup>®</sup> laser therapeutic device with KL-ILBI-525-2 and KL-ILBI-365-2 laser emitting heads for intravenous laser blood illumination in addition to kit to implement the most effective version of the combined technique ILBI-525 + LUVBI<sup>®</sup>.

Our specialists provide consultations via email concerning the application of low-level laser therapy in veterinary medicine: [7652612@mail.ru](mailto:7652612@mail.ru)